

Instruction manual



- 1.Name: derma pen
- 2.Material: ABS
- 3.Rated voltage:5V
- 4.Rated current:1A
- 5.Machine color:White, gold, silver, pink
- 6.Press and hold the start button for 3-5 seconds to turn on, and press and hold for 3-5 seconds to turn off.
- 7.When starting the machine, click on upshift (a total of 9 gears), double-click to turn on and off lighting.
- 8.“ON”Click: Work once, long press to continue working.

Features: Portable and Rechargeable

When charging the machine, please do not start the machine to use it.

How to use



Choose the right accessories



Long press for two seconds to open



Adjust the intensity



Working

Product include:

1* derma pen

1*USB charging cable

1* manual

1* box