

# Instruction manual



## Feature

1.7 Colors led Face Mask.

2. When in use to fit the face and feel comfortable.Use it,you can sit or lie down.

3.high quality non-heat producing LED lights.You can adjust its brightness by using the control.

4.Working time function: The working time of the led Face Mask is 5~60 minutes, you can adjust

the time you want to use, and enjoy your working time easily and comfortably.

( Generally we recommend Use about 15~30 minutes ) .

# Instructions



Wash your face and wait till it dries



Connect the USB cable and adapter connect LED mask and remote control



Press the "Power Switch" button and then set the time, mode, color and intensity



Put on your mask, Keep your LED mask close to your facial skin

Product Name: Led face mask

Material: ABS

Voltage: 100-240V, 50/60Hz

Colors of Lights: 7 Colors